



BECTON DICKINSON (THAILAND) CO.,LTD
24TH FL., THE OFFICES AT CENTRAL WORLD 999/9
RAMA I
PATHUMWAN
BANGKOK
THAILAND
10330
TEL +662 646 1800

Service Report

TO:	FOR:
Kanchana Sonjai Phrachomklao Phetchaburi Hospital 53 Rodfai Rd., T. Klongkrachang Phetchaburi, 76000 032-425505	Work Order Number#: WO-04237279 Case Number#: 02934307

WORK ORDER DETAILS			
Top-Level		Order Type	Field Service
Serial Number	E34297300589	Billing Type	Time and Material
Product	342973 - FACSCalibur Cytometer 3 Color Basic IVD	Customer Down	No
Work Order Task Code	T113 Repair / Troubleshooting Visit		
SERVICE REQUESTED / OBSERVED SYMPTOM / ERROR CODE, BY THE CUSTOMER			
Pressure switch can't turn on calibur			
IDENTIFIED PROBLEM			
Pressure switch can't turn on calibur			
WORK PERFORMED			
Replace new pressure on off completed, check pressure in and out pass , check sample pressure and sheath pressure pass m run beads location check pass , the instrument working normally			
Replace part , the instrument working normally			

PARTS USED

The customer acknowledges that the services described above have been provided to customer's satisfaction. Unless otherwise agreed in a valid agreement between the customer and BD, the customer agrees to pay BD's invoice for the services described above at BD's standard rates and subject to BD's standard terms and conditions.